

 **EM 650**  
**ELECTROCARDIOGRAPH**  
**6 Channels**





### TECHNICAL INFORMATION

- 6-channel ECG and adjustable for 2,3 and 4 channels
- Simultaneous displaying of 12 leads ECG waveform
- User friendly
- Color display and 5" touch screen
- Ability to record up to 6 channels simultaneously with the paper size of 110mm
- Data storage of 100 records and capability to copy of the last stored data (Save & Copy)
- Rechargeable lithium battery with 8 hours of operation and more than 150 record
- Filter types: EMG (reject muscular noise), Drift (reduce signal oscillations), HUM (50/60 Hz), Low Pass filter(25, 35, 75 Hz)
- Upgradeable software
- Comply with IEC 60601-2-25
- 16-key keyboard
- Protected against electroshock
- Displaying Rhythm-lead waveform on the screen
- Lightweight and compact design
- Different modes of recording (Manual/ Auto/ Rhythm)



**EM 650**  
**ELECTROCARDIOGRAPH**  
**6 Channels**

